


from Solutions in Critical Care

Ventilator Circuit

72" Disposable Ventilator Circuit with Smooth Bore Tubing

A high-quality, low-cost ventilation circuit compatible with adult and pediatric Bi-Level and CPAP systems.



Technical Specifications

- 1 each 180 cm smooth bore tubing - ID:19.5mm(30.5mm) OD:23.5mm(30.5mm)
- 2 each PVC straight connector - 19M-22F
- 1 each Straight Connector with tethered Sampling Port - 22M/15F-22M/15M
- 1 each 2m sampling line - ID:3.5mm OD:5.3mm; Wall Thickness:1.8mm
- 1 each hose hanger (certain configurations only)

Item ID	Description	U/M
TST600	72" Ventilator Circuit with Smooth Bore Tubing	20/Bx

Contact Your Local Distributor



Solutions in Critical Care

40 Boroline Road, Allendale, New Jersey 07401
www.solutionsincriticalcare.com • (201) 512-1100



Contracted Supplier

PP-DS-743