



DiMax & DiMax ZERO

NIV TOTAL-FACE MASKS

Transformational Masks to Meet Patients' Needs

- For use in mechanical ventilation in PSV, BiLevel-BIPAP, CPAP
- Innovative design features separate inspiratory and expiratory connections to reduce CO2 rebreathing (*DiMax ZERO only*)
- Lightweight, flexible masks feature a 0.3 mil silicone membrane allowing them to comfortably fit a wide range of face shapes and sizes without causing skin injuries.
- Five-point headgear ensures the mask stays in place, optimizes patient/ventilator synchronization and reduces leakage.
- Included connectors allow the mask to be transformed into a vented or non-vented mask or used with nasogastric tubes and fiberoptic bronchoscopic probes without leakage or ventilation interruption. Suitable for single and double-limb circuits.

Multiple Included Attachments

- Five-point headgear
- Non-Vented 22F Connector
- Non-Vented 22F Connector with FBS Port
- Vented 22F Connector with AAV
- Nasogastric (NGT) Port
- 22F - 22F Adapter (*DiMax ZERO only*)
- 22M - 22M Adapter (*DiMax ZERO only*)
- 22F Cap (*DiMax ZERO only*)



DiMax ZERO



DiMax



DiMax

Item ID	Description	U/M
700/7309US	DIMAX Total Face Mask, Small, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port)	5/box
700/7310US	DIMAX Total Face Mask, Medium, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port)	5/box
700/7311US	DIMAX Total Face Mask, Large, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port)	5/box
700/7312US	DIMAX Total Face Mask, XL, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port)	5/box
700/7313US	DIMAX Total Face Mask, XXL, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port)	5/box



DiMax ZERO

Item ID	Description	U/M
700/7315US	DIMAX ZERO Total Face Mask, Small, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port) two adapters (22-22M, 22-22F) and a 22F cap	5/box
700/7316US	DIMAX ZERO Total Face Mask, Medium, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port) two adapters (22-22M, 22-22F) and a 22F cap	5/box
700/7317US	DIMAX ZERO Total Face Mask, Large, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port) two adapters (22-22M, 22-22F) and a 22F cap	5/box
700/7318US	DIMAX ZERO Total Face Mask, XL, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port) two adapters (22-22M, 22-22F) and a 22F cap	5/box
700/7319US	DIMAX ZERO Total Face Mask, XXL, headgear, ngt port, three elbows (non-vented 22F, vented 22F with AAV, non-vented 22F with FBS port) two adapters (22-22M, 22-22F) and a 22F cap	5/box

