




Solutions in Critical Care

CPAP System

WITH INSPIRATORY/EXPIRATORY HEPA FILTER


Provide Maximum Patient Benefits & Infection Control with One Device

The Solutions in Critical Care (SCC) CPAP System is designed to deliver optimal FiO_2 levels while conserving oxygen resources. Known as the disposable CPAP with the lowest oxygen consumption, the SCC CPAP is also the only disposable CPAP with an inspiratory/expiratory filter for infection control. With multi-function access ports, the SCC CPAP is cost-efficient, easy-to-use, delivers maximum patient results and minimizes the risk of infection for patients and healthcare providers.



SCC CPAP Features

- Lowest oxygen consumption of any disposable CPAP: 5-10 Lpm
- Filter captures 99.997% of 1 μm particles (TB and Bacteria) and 99.993% of .28 μm particles (Viruses) that are inhaled and exhaled
- Port for oxygen tubing connection to boost FiO_2
- Port for nebulizer tubing connection to deliver medications
- Calibrated pressure manometer
- Low cost to buy and operate



www.solutionsincriticalcare.com
(866) FOR-SCC1 or (866) 367-7221