



Solutions in Critical Care



Panther 3

PERFORMANCE-LEVEL INVASIVE ICU VENTILATOR

3-IN-1 ICU Ventilator

The Panther 3 is a performance-level ICU ventilator offering high-performance invasive and non-invasive ventilation as well as high-flow O₂ therapy as a standard feature.

The Panther 3 ventilator can support neonates to adults and is designed to be operational in the toughest environments. With a 12.1" screen, the P3 is smaller and lighter making it a more portable device versus the Panther 5. As with all Panther ventilators, the Panther 3 does not require compressed gas or proprietary connections and was designed to require minimal maintenance.



Panther 3 Invasive Ventilator Highlights

- User-friendly 12.1" touchscreen with large, clear numbers and intuitive menus
- Neonates to adults
- Independent from compressed air
- Minimal maintenance
- Optional features available

- Invasive Modes: AC, SIMV, SMART, SPONT, APRV/BiLevel
- NIV Modes: AC, SIMV, SMART, SPONT, NIPPV
- 200 LPM Leak Compensation; 350 LPM Max
- NIV & HFOT Capabilities
- Integrated Oximetry, Capnography & Nebulizer

**LEARN
MORE**

Access the
Datasheet



Solutions in Critical Care

866-367-7221

www.solutionsincriticalcare.com

Solutions in Critical Care is the Exclusive U.S. Distributor of Origin Medical Panther Ventilators

MRK-SCC-P3-25-2.01