



inVent
RESPIRATORY

HFT150 Heated High-Flow System

Feel the Difference ... in Every Breath

Making it easier for clinicians to make a difference

The HFT150® Heated High-Flow System is a best-in-class high-flow therapy solution designed as a sleek, fully integrated device. This all-in-one system delivers a precise mix of warmed and humidified oxygen and air at flow rates ranging from 2 to 60 LPM for spontaneously breathing patients.



ALL-IN-ONE SOLUTION

- Simple, user-friendly interface
- Integrated air/O₂ blender
- Integrated heated humidifier
- Built-in nurse call alarm
- Built-in battery
- Built-in SpO₂ monitoring via Pulse Oximeter accessory
- Optional cart increases mobility and convenience



CLINICAL EXCELLENCE

- Independent temperature, flow and FiO₂ setting
- Direct FiO₂ adjustment
- Active humidity & temp control
- Servo Auto Flow (SAF) delivers precise set flow regardless of leak or back pressure
- Comprehensive alarms for patient safety



TIME-COST SAVINGS

- Simplified cleaning process
- Quick turnaround time for high-asset utilization
- Cost-effective replacement kits
- Low maintenance costs

**Contact Us
for A Demo**



HFT150 Heated High-Flow System



Item ID	Description	U/M
300-90001	HFT150 Hospital System (Power Adapter, Internal Battery and Single Disposable Kit)	1/BX
300-90009	Inlet Filter	10/BX
300-90010	Power Adapter	1/BX
300-90011	Pulse Oximeter with Custom Connector for HFT150 Interface	1/BX
300-90012	Pulse Oximetry Sensor – Small	1/BX
300-90013	Pulse Oximetry Sensor – Medium	1/BX
300-90014	Pulse Oximetry Sensor – Large	1/BX
300-90017	HFT Cart	1/BX
300-90032	HFT150 Disposables Kit (Heated Patient Circuit, Tubing Adapter, 22mm Adapter, and Humidifier Tub)	10/BX