



Solutions in Critical Care



# Panther 5N

Dedicated  
NIV Ventilator

## THE FUTURE OF HOSPITAL NIV IS HERE

### The New Gold Standard in NIV

The Panther 5N is a high-performing, reliable solution that takes NIV success to a new level by delivering

- Industry-Leading Synchrony
- Unmatched Leak Compensation
- Superior Patient Comfort

The ventilator's proprietary algorithm, state-of-the-art proportional flow valve and powerful built-in dual-stage blower can handle any leak. The P5N quickly and automatically responds to meet the patient's demand, reducing patient effort and enhancing comfort.

### Sought-After Features & Functionality

- Active Standby and AVAPS
- NIPPV and nCPAP for neonates
- User-friendly 15" touchscreen with large, clear numbers and intuitive menus
- Use with a variety of circuits and masks, even vented masks
- Dedicated NIV with HFOT and invasive capabilities
- Low cost of ownership and minimal maintenance
- Neonates to adults

\* Indicates Modes with Ramp



Features  
**ACTIVE  
STANDBY**

- NIV Modes: PCV\*, S/T\*, AVAPS, SPONT, CPAP
- 200 LPM Leak Compensation
- Up to 340 LPM Max Flow
- HFOT & Invasive Capabilities
- Integrated Oximetry, Capnography & Nebulizer
- Invasive Modes: A/C-PC, S/T, A/C-PRVC

**LEARN  
MORE**

Access the  
Datasheet



Solutions in Critical Care

866-367-7221

[www.solutionsincriticalcare.com](http://www.solutionsincriticalcare.com)

Solutions in Critical Care is the Exclusive U.S. Distributor of Origin Medical Panther Ventilators