



EAGLE™

Disposable NIV Full-Face Mask

Efficient, Affordable Solution for Short-Term NIV Therapy

The Eagle 1 NIV Full-Face Masks are a low-cost solution to providing short-term NIV therapies. The disposable mask configurations are non-vented and can be used with any ventilator circuit with an exhalation valve. Use with transport ventilator circuits and short-term non-invasive mechanical ventilation.

- Low dead space
- Lightweight, easily attached headgear
- Three color-coded sizes
- Water-clear mask dome
- Integrated chin cup
- 22mm female port



Item ID	Description	Mask Dead Space Volume	U/M
110-115019	Eagle 1 Mask, HG, Swivel Port, Large (Blue)	145 ml	8/Box
110-115020	Eagle 1 Mask, HG, Swivel Port, Medium (Green)	110 ml	8/Box
110-115021	Eagle 1 Mask, HG, Swivel Port, Small (Yellow)	80 ml	8/Box